


N95 FLAT FOLD PARTICULATE RESPIRATOR







FOLDED | MOLDED | FOLDED

FastFit combines the convenience of a flat fold with the comfort of a contour molded mask. FastFit opens quickly and easily and contours comfortably to your face. It then simply refolds for storage.

PACKAGED FOR CONVENIENCE

Individually wrapped for hygiene, it'll be handy and clean when you need it. Packed in a dispenser style box for quick and easy access when and where you need them.

EASY FIT, MORE COMFORT

Contour molded nose bridge with a soft foam nose cushion provides the wearer an easy comfortable fit. Softspun[®] lining gives added comfort and durability. There's no metal noseband to adjust and the seal is maintained without pressure points.

DURA-MESH® SHELL LASTS LONGER

The exclusive Dura-Mesh shell protects the filter media so it stays cleaner-looking longer. The shell also resists collapsing in hot, humid environments so the respirator holds its shape longer. The result is fewer respirators used and lower costs.

FEATURES

- Flat fold convenience with contour molded comfort
- · Quick and easy to put on, unlike any other flat fold
- Individually packaged for hygiene
- Packed in a dispenser box for easy access
- Proven durability and cost savings
- Dura-Mesh shell resists collapsing in heat and humidity.
- Facepiece is naturally contoured without pressure points.
- Straps are securely attached so they don't break at the seal
- Softspun® lining for increased comfort and durability.
- Soft foam nose cushion for added comfort.
- Meets heat and flame resistance in accordance with ANSI/ISEA 110-2009 Section 7.11.1.
- 100% PVC Free for a greener alternative.



N95 APPLICATIONS:



N95 FLAT FOLD PARTICULATE RESPIRATOR

| Product | Size | Valve | Strap | Display | Case |
|----------|------|-------|----------|----------|----------|
| Code | | Type | Material | Quantity | Quantity |
| 2112 N95 | M/L | None | Latex | 20/Box | 10 Boxes |

WARNING TO USER

- 1. Follow all instructions and warnings on the use of this respirator and wear during all times of exposure. Failure to do so will reduce respirator effectiveness, wearer protection and may result in sickness or death.
- 2. Before use, the user must first be trained by the employer in proper respirator use, in accordance with applicable safety and health standards, for the contaminant and exposure level in the assigned work area.
- 3. The particulates and other contaminants which can be dangerous to your health include those which you cannot see, taste, or smell,
- 4. This product has not been sold with warnings or use instructions for personnel involved in healthcare or related situations, where there may be the possibility of contact with disease or biological hazards. If you are considering such uses, first call the Moldex Technical Service Department: +1 (800) 421-0668 or +1 (310) 837-6500 ext. 512/550.

LIMITED WARRANTY IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Moldex will be met by products sold. Exclusive Remedies: damages for the breach of this limited warranty are limited to the replacement of such quantity of Moldex, products proved to be defectively manufactured. Except as provided above, Moldex shall not be liable or responsible for any loss, damage, or liability, direct, indirect, incidental, special, or consequential, arising out of sale, use, or misuse, or the inability to use products by the user.

MOLDEX TECHNICAL ASSISTANCE

For more information call +1 (800) 421-0668 or +1 (310) 837-6500 ext. 512/550 or E-mail: sales@moldex.com, tech@moldex.com or visit www.moldex.com.



Moldex-Metric, Inc. CULVER CITY, CALIFORNIA, U.S.A. +1 (800) 421-0668



HESE RESPIRATORS ARE APPROVED ONLY IN THE FOLLOWING CONFIGURATIONS

| | TC- | PROTECTION1 | RESPIRATOR | CAUTIONS AND LIMITATIONS ² | | | | |
|--|----------|-------------|------------|---------------------------------------|--|--|--|--|
| | | | 2112 N95 | | | | | |
| | 84A-5223 | N95 | Х | ABCJMNOP | | | | |

1 Protection

N95 Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply

- 2 Cautions and Limitations

- 2. Catutoris and continuations

 A- Not for use in atmospheres containing less than 19.5 percent oxygen.

 B- Not for use in atmospheres immediately dangerous to life or health.

 C- Do not exceed maximum use concentrations established by regulatory standards.

 J- Failure to properly use and maintain this product could result in injury or death.

 M- All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.

 N- Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configurations as specified by the
- Reversi equations, and for maintenance manuals for information on use and maintenance of these respirators.
 Paler to user's instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
 MIOSH does not evaluate respirators as surgical masks.

USE AGAINST:

1. Particulate aerosols free of oil.

DO NOT USE AGAINST:

- 1. Concentrations of contaminants which are unknown, or are immediately dangerous to life or health.
- 2. Concentrations of particulates which exceed the maximum use concentration or 10 times the OSHA Permissible Exposure Limit, whichever is lower.
- 3. Gases, vapors, asbestos, paint spray, sandblasting or particulate materials which generate harmful vapors.
- 4. Oil-based mists

RESTRICTIONS

- 1. This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
- 2. Do not remain in contaminated area if any physical distress occurs, for example breathing difficulty, dizziness or nausea.
- 3. Leave contaminated area and replace respirator if it is damaged, distorted, a proper fit can not be obtained or breathing becomes difficult.
- 4. Prior to each use, carefully inspect the entire respirator, including filter media and strap/weld attachment area for tears and damage. Staple or strap/weld attachment perforations do not affect respirator performance.
- 5. Do not alter, modify, or abuse this respirator.
- 6. Store unused respirators in box/bag in a clean, dry, non-contaminated area.
- 7. Dispose of respirator according to your employer's policy and local regulations.
- 8. Dispose no later than thirty days after first use.
- 9. Use respirator before the "use by" expiration date printed on box/bag.
- 10. If used for welding, wear appropriate eye and face protection.

DISTRIBUTED BY:



MOLDEX-METRIC, INC.

10111 Jefferson Blvd. Culver City, CA 90232

TEL: +1 (800) 421-0668 or +1 (310) 837-6500

FAX: +1 (310) 837-9563 F-mail: sales@moldex.com www.moldex.com



Tel: +1 (800) 421-0668, Ext. 550/512 Fax: +1 (310) 837-9563





