



Name : 6-Head Stick

Dimension : 8cm x 8cm x 2cm

Material : Red - ABS Plastic
Black – Silicone
(Anti-germ)


Massage : Hard-Percussion

Applicable Model : CM-05 & CM-07

Characteristic

- Wide area massage
- Where strong massage is needed
- Where there is more flesh than bones


Application



Shoulder



Back



Waist



Feet



Calf



Name : Point Stick

Dimension : 3.3cm x 3.3cm x 4.2cm

Material : ABS Plastic

Massage : Hard-Percussion

Applicable Model : CM-05 & CM-07

Characteristic

- Focused massage
- Acupuncture-like effect
- Helps blood circulation

Application



Applied Image



Applied Image



7.1 cm

7.1 cm

Name : Air-Cushion Stick

Dimension : 7.1cm x 7.1cm x 4.2cm

Material : Red - ABS Plastic
Black - Silicone





Massage : Soft-Percussion

Applicable Model : CM-05 & CM-07

Characteristic

- Air-pumping type massage stick
- Where there is more bone than flesh
- Where bone is near like wrists, feet, fingers, knees

Application





4.2 cm

4.2 cm

Name : Facial Massage Stick

Dimension : 3.3cm x 3.3cm x 4.2cm

Material : Red – Silicone
Silver – Stainless Steel

Massage : Vibrating

Applicable Model : CM-07


Characteristic

- . Facial massage stick
- . Stainless-steel applied

Application



Applied Image



Applied Image



5.2 cm



5.2 cm

Name : Body Oil Massage Stick

Dimension : 5.2cm x 5.2cm x 1.6cm

Material : Silicone (Anti-germ)

Massage : Soft-Percussion
Vibration

Applicable Model : CM-05 & CM-07

Characteristic

- Works best with body gel, crème, or body oil
- Compatible between Head and Grip attachment
(Head attachment for percussion massage)
(Grip attachment for vibration massage)


Application



1. Rub body gel



2. Attach the stick



3. Massage



5.2 cm

5.2 cm

Name : Scalp Massage Stick

Dimension : 5.2cm x 5.2cm x 1.6cm

Material : Silicone (Anti-germ)

Massage : Soft-Percussion
Vibration

Applicable Model : CM-05 & CM-07

Characteristic

- Scalp massage stick
- Compatible between Head and Grip attachment
(Head attachment for percussion massage)
(Grip attachment for vibration massage)

Application



Applied Image



Attached to Head
-. Percussion massage



Attached to Grip
-. Vibration massage